

THE T-SPOT.TB BORDERLINE RESULT

The vast majority of T-SPOT.TB test results are either Positive or Negative; however, a small percentage of test results can be Borderline (equivocal), where the higher of (Panel A minus nil control) and (Panel B minus nil control) is 5, 6 or 7 spots.¹

Why does the T-SPOT.TB test have a Borderline category?

- Interferon-Gamma Release Assays (IGRAs) including the T-SPOT.TB test, are designed to measure an adaptive immune response to TB and are therefore subject to biologic variability²
- To address test variation and uncertainty of results near a dichotomous cut point, the T-SPOT.TB test has a Borderline category, which is intended to minimize false-positive and false-negative results around the test cutoff² (See Figure 1)
- King et al. reported that the use of the FDA-approved Borderline category for the T-SPOT.TB test resulted in substantially lower observed conversion and reversion rates than without its use, suggesting that the Borderline category is useful in decreasing the number of discrepant results²
- According to the Centers of Disease Control and Prevention (CDC), "Incorporation of a borderline category for the T-SPOT [.TB] test as approved by FDA increases test accuracy by classifying results near the cut point (at which small variations might affect the interpretation) as neither positive nor negative"³

What should I do with a Borderline T-SPOT.TB test result?

- When a Borderline result is obtained, retesting by collecting another sample is recommended¹
 - In a large multi-center study by King et al., the majority of Borderline results, 79.8%, resolved as clearly Positive or Negative upon retesting²
 - This study also demonstrated that 23% of Borderline results were Positive on retesting, concluding that the Borderline category is clinically relevant and useful in maintaining test sensitivity²
- Upon retesting, if the test result remains Borderline, other diagnostic tests and/or epidemiologic information should be used to help determine the TB infection status of the patient¹

How long should I wait to retest a patient after receiving an initial Borderline test result?

- There is no set guideline established for the time interval between an initial Borderline test result and a retest
- The ordering provider, having access to the patient's medical history, is best equipped to determine the appropriate time interval between tests




Figure 1. Hypothetical graph of test results from true-positive and true-negative subjects

Most quantitative tests do not perfectly discriminate between subjects with and without a given disease.⁴ To address test variation and uncertainty of results near a dichotomous cut point, the T-SPOT.TB test has a Borderline category, which is intended to minimize false-positive and false-negative results around the test cut-off.

THE T-SPOT.TB INVALID RESULT

A small percentage of all T-SPOT.TB results will not be considered clinically interpretable and will be reported as Invalid.

What is an Invalid T-SPOT.TB result?

- A T-SPOT.TB result is considered Invalid if the positive and/or nil (negative) control does not perform as expected:¹
 - The nil control is designed to control for non-specific T-cell reactivity.¹ A nil control spot count in excess of 10 spots should be considered as 'Invalid'
 - The positive control serves as an indicator of patient cell functionality.² If the positive control spot count is < 20 spots, the result should be considered Invalid, unless either Panel A or Panel B are Positive or Borderline
- Invalid results are uncommon (<1%⁵) and may be related to the immune status of the individual being tested. They may also be related to a number of technical factors, potentially resulting in high background, low mitogen, and high nil results, including:^{1,2}
 - Inappropriate blood storage conditions
 - Delay in sample transport
 - Patient specific conditions
 - Laboratory error

What should I do with a Invalid T-SPOT.TB test result?

- Invalid results are not clinically interpretable and retesting by collecting another sample is recommended²
- Upon retesting, if the test result remains Invalid, other diagnostic tests and/or epidemiologic information should be used to help determine the TB infection status of the patient²

How does a T-SPOT.TB Invalid result differ from a T-SPOT.TB Borderline result?

- Invalid results cannot be interpreted due to the failure of the test's positive and/or nil (negative) control¹
- A Borderline result is a valid test where the spot count is too close to the established test cut-off to classify the result as Positive or Negative

REFERENCES:

1. Oxford Immunotec. T-SPOT.TB Package Insert PI-TB-US-0001 V7. Abingdon, UK. February 2019
2. King TC, Upfal M, Gottlieb A, Adamo P, Bernacki E, Kadleck CP, Jones JG, Humphrey-Carothers F, Rielly AF, Drewry P, Murray K, DeWitt M, Matsubara J, O'Dea L, Balsler J, Wrighton-Smith P. T-SPOT®.TB Interferon-Gamma Release Test (IGRA) Performance in Healthcare Worker Screening at 19 US Hospitals. American Journal of Respiratory and Critical Care Medicine, 2015, Vol.192: 367-373.
3. Mazurek GH, Jereb J, Vernon A, LoBue P, Goldberg S, Castro K, Committee IE, Centers for Disease Control and Prevention. Updated guidelines for using Interferon Gamma Release Tests to detect *Mycobacterium tuberculosis* infection - United States, 2010. MMWR Recomm Rep 2010; 59: 1-25.
4. Coste J, Pouchot J. A grey zone for quantitative diagnostic and screening tests. *Int J Epidemiol.* 2003; 32(2):304-313.
5. Rego K, Pereira K, MacDougall J, Cruikshank W. Utility of the T-SPOT®.TB test's borderline category to increase test resolution for results around the cut-off point. *Tuberculosis.* 2018;108:178-185.doi:10.1016/j.tube.2017.12.005.

T-SPOT®.TB

 **Oxford
Immunotec**